

Test system SIMICON EN-DUO-PROT

for endoscope washer disinfectors type MEDIVATORS (single chamber / dual chamber)

Area of application:

The test system SIMICON EN-DUO-PROT is designed for the evaluation of protein removal in WD for flexible endoscopes. The test system contains the indicator SIMICON RI in endoscope process challenge devices simulating endoscope channels. The indicator SIMICON RI in the SIMICON PCD is in terms of its resistance comparable to the test piece in the German DGKH Guideline as well as compliant to ISO/TS 15883-5.

Testing procedure:

Proceed likewise for each chamber.

1

Placement of the endoscope PCDs

Take the PCDs out of the aluminium pouch and connect them to the respective rinsing nozzles.

Note: The PCDs already contain the indicator SIMICON RI.
Please do not open the receptacles nor unroll the PCDs.

2

Starting the Process

Choose the disinfection programme and start the washer disinfectant.

IMPORTANT: The detection of residual protein calls for the abortion of the process after the cleaning cycle.

3

Test report

Enter all requested data in the test report. Please provide your address, telephone and fax number.

Note: Mark the tested chamber on the report.

4

Removing the indicators

Note: For the following steps use one-way gloves.

Disconnect the PCDs and place them in the self-adhesive tyvec-pouches provided. Label the pouches with the corresponding chamber and PCD numbers (s. drawing).

Shipping

Return the test system to SIMICON. Please check first for completeness.

As soon as the assessment is completed, SIMICON will return the test report including test results to you for your documentation.

