



DOUBTS?

Play it safe!

**SIMICON RI, the optimized system
for monitoring the cleaning efficacy
of washer-disinfectors.**

**Easy handling.
Clear results.**



Unique precision!

Special PCDs for simulating standard instruments,
narrow gaps and MIS-instruments.

SIMICON 

Exceptional Hygienic Safety




WITHOUT A DOUBT!

SIMICON RI is the **reliable cleaning indicator** for validation and routine monitoring of cleaning and disinfection processes. Applicable for MIS and standard instruments as well as for endoscope and ultrasonic cleaning. Special PCDs simulate distinguished requirements of the cleaning process. **If the operating parameters are set correctly, the test soil will be removed completely. Failures will be identified reliably.** SIMICON RI for fast and clear visual assessment of the cleaning result.

Advantages at a glance:

- Easy handling
- Clear evaluation of the cleaning efficacy
- Optimized simulation of standard instruments, narrow gaps and MIS-instruments
- Conformity according to ISO 15883
- PCDs for various systems

Certainty at first glance:



Before the cleaning process:
Indicator with test soil.



After the cleaning process:
Test soil completely removed.
Easy assessment.



SIMICON 

Exceptional Hygienic Safety



NOW

Order product information and prices:

By mail: info@simicon.de

By Fax: +49 (0) 89/673366-22

Company / Hospital / Practice	
Department	
Forename / Surname	
Road / House Number	
Postcode	Town
Phone / Fax / E-Mail	

SIMICON GMBH · Sigmund-Riefler-Bogen 19 · 81829 München
Phone: +49 (0) 89 / 67 33 66-44 · www.simicon.de

SIMICON 

Exceptional Hygienic Safety